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Contents of these Operating Instructions

1. Models
2. Technical data / Drawings
3. Intended uses
4. General safety information / Warranty
5. Transport
6. Installation
7. Interior fittings
8. Storage amounts / Storing substances together
9. Ventilation
10. Grounding
11. Markings
12. Decommissioning after a fire
13. Care and maintenance
14. Disposal

1. Models

- 1.1 Type 90 – safety cabinets as per EN 14470-1
- 1.2 Type 30 – safety cabinets as per EN 14470-1
- 1.3 G 30 – gas cylinder safety cabinets as per EN 14470-2

2. Technical data / Drawings

see Annex 1

3. Intended uses

- **Type 30 / Type 90 safety cabinets are suitable for the storage of flammable liquids at the workplace while observing the relevant national rules and regulations.**
- **G 30 gas cylinder safety cabinets are suitable for the storage of gas cylinders at the workplace while observing the relevant national rules and regulations.**

4. General safety information / Warranty

- For the operation of your safety cabinet you must observe the national rules, regulations and laws that are applicable.
- Before starting to use your safety cabinet you must carry out a functional check in order to recognise any damage that may have been caused by transport, importation or improper installation!
- Safety devices must be checked at regular intervals – this should only be carried out by our trained specialists! This guarantees proper and safe long-term storage
- Warranty: in accordance with our General Terms and Conditions. See also item 13.
- Please also observe your duty of exercising due care – any damage must be remedied and spilt liquids removed immediately!
- The area in which the doors open must be kept free and clean.
- Doors and drawers should only be open when items are being removed from or placed in the safety cabinet – keep the doors/drawers closed!
- Please observe the maximum amounts that can be stored, shelf loading, etc., as per Annex 1 and the markings on the door/drawer.
- The safety cabinets are not intended for the storage of aggressive substances (e.g. acids and alkalis, corrosive gases) as these could corrode the interior fittings and the exhaust air extraction system! Our special acid and alkali safety cabinets should be used for this.
- Before storing substances (e.g. reagents) you must check that the cabinet surfaces are resistant to them.
- To comply with EN 14470-1 the collection trough must be able to contain a minimum volume of 10% of the contents of the containers stored above it, or at least 110% of the contents of the largest single container – do not store too large containers!
- See 7.4 perforated cover plates

5. Transport

storeLAB® safety cabinets must only be lifted from beneath by using a suitable means of transport; they must be transported upright. Tilting during lifting and transport must be avoided at all costs. Tilting could cause damage that could limit or destroy the functioning of the built-in safety devices. You must also avoid jolting the cabinet when unloading it.

Always lift the cabinet from the side or the rear and never from the door side. The carrying surface of the hand truck, fork-lift truck or similar must have at least the same area as the base of the safety cabinet.

6. Installation

The safety cabinet must be installed, operated and maintained in such a way that employees are safeguarded from the risks of fires and explosions.

See assembly instructions (inside cabinet door).

7. Interior fittings

7.1 Shelves

- Installed according to assembly instructions (inside cabinet door)

7.2 Pull-out shelves

- Already built into cabinets

7.3 Shelves

- Already built into cabinets

Please observe the maximum shelf load (see Annex 1)!

If the shelves are not installed in accordance with the assembly instructions then the maximum shelf load cannot be guaranteed.

7.4 Collection trough / Perforated cover plate

The collection troughs in the base of the storeLAB® safety cabinet or the pull-out shelves are intended to collect any liquid leakages inside the cabinet and are not additional support surfaces for stored containers (see EN 14470-1).

If a perforated cover plate is used for the collection trough / pull-out shelf then this can be used as additional support surfaces for stored containers.

8. Storage amounts / Storing substances together

Valid national regulations must be observed.

9. Ventilation

All safety cabinets have an air inlet and an air outlet connection. The outlet connection (suction side, connection diameters attached annex 1) is connected to the exhaust system. The exhaust air is led away through ducting to discharge into the atmosphere at a non-hazardous location. Inlet air is taken from inside the room. Please observe the valid national rules, regulations and laws concerning the stipulated rate of air exchange and exhaust air removal.

Safety cabinets as per EN 14470-1 can also be operated without connection to a ventilation system. (In Germany please observe the requirements of TRbF 20-L; elsewhere the relevant national rules, regulations and laws.)

Please remember the following points for operation without connection to a ventilation system:

- The cabinet must be grounded to avoid electrostatic charges.
- The area 2.5 m around the safety cabinet up to a height of 0.5 m above floor level is to be regarded as being a Zone 2 explosion- endangered area.
- In accordance with the relevant accident prevention regulations (BGV A8), this Zone 2 explosion- endangered area must be clearly and visibly marked with the following hazard symbol:

The locations of the air inlet and outlet connections are shown under Annex 1 (see drawings).

10. Grounding

In accordance with BGR 132 (D), all safety cabinets are equipped with a grounding connection on the cabinet roof as standard. If liquids are transferred within the cabinet then the interior of the cabinet must be additionally grounded.

11. Markings

Markings as per EN 14470-1 / 14470-2 / TRbF 20, Annex L / ISO 3864.

12. Decommissioning after a fire

Depending on the duration of the fire and the fire load, an explosive vapour-air mixture can be formed inside any safety cabinet. Before opening the storeLAB® safety cabinet all sources of ignition, naked flames, etc., must be removed from a distance of 10 m around the safety cabinet.

- Inform the fire department about the contents of the cabinet and the period mentioned below.
- Only open the cabinet approx. 24 h after the fire has been extinguished, or wait for a period that is at least six times as long as the duration of the fire.
- Check that the outer surface of the cabinet is cold. If it is still hand-warm then the waiting period must be extended.
- Avoid under all circumstances the use of cutting tools that produce sparks and burning tools.
- Ensure that suitable fire-fighting agents are at hand during the opening procedure.
- Wear suitable protective clothing!

13. Care and maintenance

The storeLAB® safety cabinet is a technical safety device according to the relevant German and European safety regulations (complies with EEC directive). To ensure its protective function in the event of a fire it must be checked at regular intervals. The German safety regulations stipulate that the safety devices must be checked by a specialist at regular intervals. The interval between safety checks is currently one year.

13.1 Daily checks

- Escaped liquids must be removed immediately and disposed of properly.
- Check the storeLAB® safety cabinet regularly for mechanical damage, dirt and corrosion.

13.2 Monthly checks

The monthly check covers the following points:

- | | |
|--------------------|---|
| 1. Door closing: | Open the door and check that it closes properly. |
| 2. Ventilation: | Test the efficiency of the ventilation system with a piece of wool, meter, etc. |
| 3. Sealing strips: | Check that the sealing strips in the door groove and the side(s) of the door(s) are correctly seated. If any damage is visible the sealing strips must be replaced immediately. |

14. Disposal

The storeLab® safety cabinet can be completely dismantled. Individual components, e.g. metal, insulation panels, door locks, etc. can be recycled separately. National and local waste disposal regulations must be observed. To preserve resources the safety cabinet components or the safety cabinet itself should't be disposed of as domestic waste.

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**storeLAB® Safety Storage-Cabinet model 1400 / model 1200 / model 900 /
model 600 / model SIS-FAS 1550 / model SIS-FAS 900
Underbench cabinets model 1400 FT / model 1400 SF / model 1400 SFSF /
model 1400 FTSF /
model 1100 FT / model 1100 SF / model 1100 SFSF / model 1100 FTSF /
model 880 FT / model 880 SF / model 880 FT SF / model 880 SF SF /
model 600 FT / model 600 SF
model 1400-500FT / model 1400-500 SF / model 1400-500 FT SF / model 1400-500SFSF
model 1100-500 FT / model 1100-500 SF / model 1100-500 FT SF /
model 1100-500 SF SF /
model 880-500 FT / model 880-500 SF / model 880-500 FT SF / model 880-500 SF SF /
model 600-500 FT / model 600-500 SF /
Safety cabinets for documents (DKS) model 1200 / model 1400
Chemical cabinets with safety box model 950-SiB 30 / model 950 GL-SiB 30 /
model 950-SiB 60 / model 950 GL-SiB 60 / model 950-2FAS-SiB 30 / model 950-2FAS
GL-SiB 30 / model 950-2FAS-SiB 60 / model 950-2FAS GL-SiB 60**

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to which this declaration relates is in conformity with the relevant provisions of the following stand-
ards or other normative documents.

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DIN EN 14470-1 (2004-07)

Fire safety storage cabinets – Part 1: Safety storage cabinets for
Flammable liquids

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DIN/ISO 3864 Safety colours and safety signs

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Applied nation-
al standards
(germany)

DIN 4844-1 Safety colours and safety signs

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Following the provisions of Directives (Germany)

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and the regulations of the Statutory Accident Insurance Institution

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The operation company is responsible for installing the products in accordance to national regulations
which may apply in respect to the installation, use and also for regular inspections.

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to which this declaration relates is in conformity with the relevant provisions of the following standards or other normative documents.

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EG-Maschinenrichtlinie 2006/42/EG, L'157/65 vom 9.6.2006
EEC-Machinery Directive 2006/42/EEC
EN 292 Sicherheit von Maschinen; Grundbegriffe, allgemeine
Gestaltungsgrundsätze
EN 292 Safety of machinery – Basic concepts, general principles
For design

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







SiS Typ 30 / 1200

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SiZ G30

Technical drawing of a rectangular ventilation duct. The overall dimensions are 1190 mm in width and 845 mm in height. The duct has a top flange with a height of 293 mm and a bottom flange with a height of 827 mm. The distance between the top and bottom flanges is 554 mm. The duct is labeled with 'Zuluft' (supply air) on the left and 'Abluft' (extract air) on the right. The duct has a diameter of $\varnothing 80$ mm. The duct is supported by a base with a width of 162.5 mm and a height of 80 mm. The duct is labeled with 'Erdschleim' (earthworm) on the right side.

[illegible]

Technical drawing of a rectangular ventilation unit. The overall dimensions are 1400 mm in width and 598 mm in height. The unit features a 'Zuluft' (supply air) inlet on the left and an 'Abluft' (extract air) outlet on the right. The 'Abluft' outlet is connected to an 'Erdungsklemme' (grounding terminal). The unit has a series of 16 circular fans along the bottom edge, with a diameter of $\varnothing 80$ mm. The 'Zuluft' inlet has a diameter of $\varnothing 25$ mm. The 'Abluft' outlet has a diameter of $\varnothing 80$ mm. The distance from the left edge to the first fan is 162 mm. The distance between the fans is 80 mm. The distance from the last fan to the right edge is 162 mm. The distance from the 'Zuluft' inlet to the 'Abluft' outlet is 293 mm. The distance from the 'Abluft' outlet to the right edge is 156 mm. The distance from the bottom edge to the 'Abluft' outlet is 428 mm.

Technical drawing of the front view of a rectangular unit. The overall dimensions are 900 mm in width and 598 mm in height. The drawing includes the following features and dimensions:

- Zuluft** (Supply Air) inlet on the left side.
- Abluft** (Exhaust Air) outlet on the right side.
- Erddungsklemme** (Earth terminal) located on the right side, with a vertical distance of 156 mm from the top edge.
- A circular feature with a diameter of $\varnothing 80$ mm is located on the right side, with a vertical distance of 428 mm from the bottom edge.
- A circular feature with a diameter of $\varnothing 25$ mm is located on the left side.
- A series of small circular features (possibly screws or mounting points) are located along the bottom edge, with a horizontal distance of 151 mm from the left edge to the first one, and 80 mm between subsequent ones.

Technical drawing of a square ventilation grille. The overall dimensions are 570 mm by 570 mm. The central circular opening has a diameter of $\varnothing 80$ mm. The distance from the center of the opening to the top edge is 196 mm. Labels indicate 'Zuluft' (supply air) on the left, 'Abluft' (exhaust air) on the top, and 'Erdungsklemme' (grounding terminal) on the right. The drawing shows a square frame with mounting points at the corners and midpoints of the sides.

Technical drawing of the front view of a rectangular enclosure. The drawing shows a rectangle with a width of 1100 and a height of 570. On the top edge, there is a circular feature with a diameter of $\varnothing 80$. Above this feature, there is a dimension of 196. Labels with arrows point to the top edge (Zuluft), the circular feature (Abluft), and a terminal (Erduungsklemme).

Technical drawing of a rectangular box. The overall width is 877 and the overall height is 570. A top flange has a width of 195. A hole in the top flange is labeled with a diameter of $\varnothing 80$. Arrows point to the top edge (labeled 'Zuluft') and the top flange (labeled 'Abluft'). A label 'Erddungsklemme' points to the right side of the box.

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







SiS Typ 30 / 1200

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SiZ G30

Technical drawing of a rectangular ventilation unit. The overall width is 1190 and the overall height is 598. The internal width is 293. The unit features a 'Zuluft' (supply air) outlet on the left and an 'Abluft' (extract air) outlet on the right. The extract air outlet has a diameter of $\varnothing 60$. The unit is equipped with 10 $\varnothing 25$ connection points along the bottom edge. The drawing also shows an 'Erdungsklemme' (grounding terminal) on the right side. The bottom edge is divided into segments with dimensions: 162.5, 80, 80, 80, 80, 80, 80, 80, 80, 80, 162.5.

Technical drawing of a rectangular ventilation duct. The overall dimensions are 1400 mm in width and 598 mm in height. The duct has a 'Zuluft' (supply air) inlet on the left and an 'Abluft' (extract air) outlet on the right. The distance from the right edge to the center of the 'Abluft' outlet is 293 mm. The distance from the top edge to the center of the 'Abluft' outlet is 156 mm. The duct is equipped with 16 circular air outlets along its length, each with a diameter of $\varnothing 25$ mm. The distance from the left edge to the center of the first outlet is 162 mm, and the distance from the center of the last outlet to the right edge is also 162 mm. The spacing between the 15 internal outlets is 80 mm. The duct is connected to an 'Erdschleife' (grounding loop) on the right side. The drawing also shows a cross-section of the duct with a diameter of $\varnothing 80$ mm.

Technical drawing of a rectangular enclosure. The overall width is 877 and the overall height is 570. A top flange has a width of 195. A circular feature on the top flange has a diameter of $\varnothing 80$. Labels indicate 'Zuluft' (supply air) on the left side, 'Abluft' (exhaust air) on the top flange, and 'Erddungsklemme' (grounding terminal) on the right side.

Technical drawing of a rectangular enclosure. The overall width is 1400 and the height is 570. The top edge features a 'Zuluft' (inlet) on the left and an 'Abluft' (outlet) on the right. The outlet is a circular duct with a diameter of $\varnothing 80$ and a height of 195. An 'Erdungsklemme' (grounding terminal) is located on the right side of the enclosure.

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LACONT Umwelttechnik GmbH